




Cryogenic Target Upgrade


Terry Tope 4.28.07



Previous Hydrogen Target Layout



Previous Hydrogen Target Inside TPC




Previous Hydrogen Target Layout



Previous Hydrogen Target Layout




Existing Hydrogen Flask & Vacuum Sphere



New Hydrogen Target Layout

- Issues
 - Flask location?
 - Minimum “throat” area for beam?
 - Low mass vacuum shell
 - Current OD > 5.5” plastic ball aperture, OK for 10” aperture
 - Change shape to fit thru ball aperture?
 - Cylinder wall 2x thick as sphere
 - Where to mount refrigerator?
 - Magnet face or plastic ball structure?
 - How to align target?
 - Possible to run flex lines so target can move with plastic ball independent of refrigerator to move out of beam
 - Large buffer tank can be mounted remotely with a hose connection
 - Do we care about vibration?
 - Safety with HV in the area



Target Upgrade


- Purchase new cryocooler
 - Current cryocooler is 20+ years old
 - New cryocooler offers improved reliability, 20k hours until first maintenance
 - Automation of manual adjustments
 - Ability to use cryocooler for N2 target using same fluid circuit as H2 by crippling cryocooler with heaters



14 W @ 20K

80 W @ 80 K

Cryomech pulse
tube cryocooler



~\$36k for total package, 14 week delivery